TRYGGHET

Design Johanna Jelinek

Design and Quality
IKEA of Sweden

<http://ikea-club.com.ua>

ENGLISH	4
DEUTSCH	8
FRANÇAIS	12
NEDERLANDS	16
DANSK	20
NORSK	24
SUOMI	28
SVENSKA	32
ČESKY	36
ESPAÑOL	40
ITALIANO	44
MAGYAR	48
POLSKI	52
LIETUVIŲ	56
PORTUGUÊS	60
ROMÂNĂ	64
SLOVENSKY	68
HRVATSKI	72

IMPORTANT! READ CAREFULLY AND RETAIN FOR FUTURE REFERENCE!

Intended usage

This first aid kit is intended for small injuries and superficial wounds such as cuts and grazes, or small burns and soft tissue injuries like sprained wrists or ankles.

When to seek further medical help

If an injury has a high risk of bacterial infection e.g. wounds from animal bites or from dirty objects. In case of internal injuries e.g. accidental poisoning or other severe injuries, large burn wounds or heavy bleeding.

If you are not sure of the severity of an injury.

This first aid kit includes:

40 plasters in 2 sizes, 10 non-alcoholic swabs, 6 skin closure strips, 1 elastic bandage (4 m) and 3 compresses (7x7 cm).

Step by step instructions

Swabs

Swabs can be used for cleansing all small wounds and can also be used for cleaning hands before treating an injury. They contain a mild and natural salt solution, but are not antibacterial.

HOW TO USE SWABS:

1. Clean the wound from the centre outwards to avoid getting bacteria into the wound.
2. Throw the swab away after use.

TIPS:

You may need several swabs depending on the size of the injury. Instead of swabs you can use soap and water or an antibacterial solution.

Plasters

The plasters in this kit are suitable for small superficial wounds. If the wound is bigger than the absorbing pad on the plaster or for gaping or deeper wounds, use closure strips or the compress and elastic bandage.

HOW TO USE PLASTERS:

1. Wash your hands with soap and water.
2. Clean the wound from the centre outwards to avoid getting bacteria into the wound. Make sure the skin is dry before applying the plaster.

TIPS:

Replace the plaster after 24 hours to avoid infection. Change the plaster if it gets wet since germs thrive in moist conditions. Small wounds heal better when exposed to air, so remove the plaster at night if possible.

If applying a plaster on the joint of the body (knee, finger etc), then apply when the joint is bent to prevent the plaster from getting stretched and falling off, or causing discomfort. Seek medical help if there is any swelling, reddening of the skin and increased pain, since it may be a sign of infection. If the skin becomes irritated it could be due to an allergic reaction, in which case you should remove the plaster. Never reuse a plaster.

Skin closure strips

Use skin closure strips to stop the bleeding of small gaping wounds. **IMPORTANT:** for deeper cuts and soiled wounds from e.g. animal bites, please seek further medical help.

HOW TO USE SKIN CLOSURE STRIPS:

1. Wash your hands with soap and water.
2. Clean the wound starting at the centre outwards to avoid getting bacteria into the wound. Let the skin dry before applying the closure strips.
3. Place the strip at one end of the wound. Close the edges of the wound as tightly as possible and secure the strip across the wound. Some wounds may require several strips, in which case you should space them 3-5mm apart.
4. We recommend you leave the strips on for up to 10 days to allow the wound to heal.

TIPS:

Try to keep the strips from getting wet. If they do get wet and come off, change them for new ones.

You may need to cut the strips to size for smaller areas, e.g. eyebrows. If so, make sure you cut them before removing the paper sheet. Seek medical help if there is any swelling,

reddening of the skin and increased pain, since it may be a sign of infection. Consult your medical advisor if you are concerned about scarring.

Compress for wounds/burns

Use for treating bleeding wounds and small burn injuries.

A: HOW TO APPLY A PRESSURE BANDAGE FOR BLEEDING WOUNDS

For wounds with heavy bleeding or if you are unsure of the severity of an injury, always seek medical help!

1. Wash your hands with soap and water.
2. Clean the wound from the centre outwards, to avoid getting bacteria in the wound.
3. Put the compress against the bleeding wound and apply pressure. The compress will not stick to the wound since it has a glossy surface.
4. Unroll some of the bandage but leave most of the roll intact. Put the roll on the compress and apply pressure.
- 5-6. Wrap the loose part of the bandage roll around the injured body part and fasten with the included hooks. Leave until the bleeding has stopped.

B: HOW TO APPLY A COVER DRESSING FOR SMALL BURN INJURIES.

IMPORTANT! If a burn injury is proportionally larger than the palm of the hand of the injured person, you should immediately seek medical help. Be extra attentive with burn injuries of children since a small injury might be proportionally large enough to require medical help.

1. Cool the wound with clean water, approx 18-20°C (64,4-68°F), for 10-15 minutes.
2. Wash your hands with soap and water.
3. Place the compress on the burn wound to protect it from dirt.
4. Keep the compress in place by wrapping the elastic bandage around it in a loose fashion.
5. Fasten with included hooks.

Elastic bandage

Use an elastic bandage to treat soft tissue injuries such as sprained or twisted ankles and knees, or extensive bruising. You can also use it to keep a compress in place (see instructions under wound/burn compress). If you suspect a fracture, please seek further medical help.

HOW TO USE AN ELASTIC BANDAGE:

1. Keep the injured body part still and elevated to minimize internal bleeding and swelling. No cleaning is needed as long as the skin is intact.
2. Wrap the bandage tightly, in the direction towards the heart of the injured person. Fasten with the included bandage hooks and leave for 15 minutes.
3. Remove the bandage for five minutes, keeping the injured part of the body raised and still.
4. Wrap the bandage tightly once again, and leave for another 15 minutes.
5. Remove the bandage.

Following the above instructions will help the healing process. You can reapply the bandage afterwards loosely yet firmly to provide some support. However it will not help the healing process.

Remember to always clean your hands before treating an injury. And don't hesitate to seek further medical help in case you are unsure of what to do in an emergency.

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